

INTERPRETATION

Health Guidance Values(BGVs) are set at a level at which there is no indication from the scientific evidence available that the substance being monitored is likely to be injuries to health. Values not greatly in excess of a HGV are unlikely to produce serious short or long term effects on health. However, regularly exceeding the HGV does indicate that exposure is not being adequately controlled. Under these circumstances employers will need to look at current work practices to see how they can be improved to reduce exposure.

For further advice contact:

Group Support Unit, Room L.2.51
Health and Safety Laboratory,
Harpur Hill, Buxton, SK17 9JN
Tel: 01298 218099
Fax: 01298 218172

Website: www.hsl.gov.uk

BIOLOGICAL MONITORING METHODS

October 2005



Method for Mercury in Urine

Hazardous Substance:

Mercury & inorganic compounds

Occupational Exposure Standard = 0.025 mg m⁻³
CAS No. 7439-97-6

Biological Monitoring Guidance value:

HGV= 20 µmol Hg/mol creatinine
Conversion: 1 µmol/mol = 1.77 µg

❑ **Sample Collection**

Time : Random, untimed urine

Equipment: Polystyrene universal container (30ml)

❑ **Description of Suggested Method**

Analysis of urine using direct nebulisation Inductively Coupled Plasma Mass Spectrometry (ICP-MS)¹. Samples should be diluted in nitric acid, and it is advisable to add gold to the samples and standards to stabilise the analytical performance. An internal standard must be added to the samples to compensate for matrix effects.

❑ **References**

1. Morton J, Mason H J, Richie K A, White M A, Comparison of hair, nails and urine for biological monitoring of low-level inorganic mercury exposure in dental workers, *Biomarkers*, 2004, 9, 1, 47-55.

2. Richardson R A, (1976) Automated method of the determination of mercury in urine. *Clin Chem*, 22, 1604-1607

3. Mason J H, Calder I. (1994) The correction of urinary mercury concentrations in untimed, random samples. *Occup. Env. Med*, 51(4),287.

❑ **Alternative Method**

Mercury in urine may also be determined using a cold-vapour mercury detector² or by cold vapour atomic absorption spectrometry³. The method consists of firstly digesting the urine sample by reduction of Hg²⁺ to Hg⁰ using stannous chloride and subsequent measurement of elementary mercury with the detector of choice.

❑ **Sample Transport to Laboratory**

At ambient temperature, samples should arrive within 48h of collection. If delay anticipated store at -20°C. Samples sent through postal system must comply with Post Office regulations.

❑ **Analytical Evaluation**

Precision

- within day <4% RSD at 180 nmol/l
- day to day <6% RSD at 180 nmol/l
- recommended precision <7.3% RSD³

Detection limit

- 5 nmol/l
(creatinine)

Calibration range

- typically 0 - 500 nmol/l

Sample Stability

- 2 days at RT, > 6 months at -20 °C

Analytical interferences

- None known

❑ **Other Information**

At moderate exposures result reflect cumulative exposure over recent weeks/ months

Half life-time

-40-60 days

❑ **Confounding Factors**

None known

Unexposed Levels

< 2 µmol/mol creatinine

Creatinine Correction

Advised, specific gravity correction less suitable

❑ **Quality Assurance**

Internal QC - must be established

External QA - TEQAS, University of Surrey (Tel: 01483 509217)